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
Nikon AX Confocal Systems

Environmental and Placement Specifications,
Power Requirements and System Component Information

05-2021

Environmental conditions AX and AXR:

Operating conditions:

Temperature: +23 ±5°C, 24hr air conditioning required.
Humidity: 70%RH or lower for AX(R), no condensation.
60%RH or lower for the used laser units and microscopes, no condensation.
Altitude: 2,000 m or lower
Pollution level: Degree 2
Overvoltage category: Category II
Indoor use only

Storage conditions:

Temperature: 0 to 35°C
Humidity: 70%RH or lower (no condensation).

Transport conditions (within 48 hours):

Temperature: -20 to +50°C
Humidity: 90%RH or lower (no condensation)

Power Supply:

AX(R) Controller:

Input: 100-240 VAC ±10%, 50/60 Hz, 6A-3A
Overcurrent protection: Fuse 250 V 10A, fast-acting, Ø5.2 x 20 mm
Electric shock protection class: Class I

LUA-S4/ S5/ S6 Laser unit:

Input 100 – 240 VAC± 10%, 50-60 Hz, 15A max.
Electric shock protection class: Class I

Microscope:

(Reference value Ti2, please match the actual value with the actual microscope).
Power Supply: 240 VAC ±10%, 50/60 Hz, < 3A

Handle the system carefully.

This system is a precision optical instrument. Handle the system with care to avoid physical shock or impact. Gently handle the cables. In particular, do not bend or pull the optic fiber cables with excessive force to prevent failure, malfunction, or breakage of fibers.


The microscope and the AX scan head are precisely adjusted. Subjecting these instruments to shock or forced relocation may result in the unevenness of acquired images.

Installation conditions.

To ensure safe and stable observations, check the following conditions when installing this system.

- To prevent the system from falling and to avoid malfunction, install each unit of this system in a flat place.
- Use the system in a clean environment. Accumulated dust or foreign particles may significantly degrade the performance.
- Use the system in the environment specified in Chapter 7, "Specifications". To prevent condensation or failure, do not use the system under conditions of high temperature and humidity.
- To prevent vibrations and subsequent degradation in image quality, install the microscope and the AX scan head on a solid surface, such as a vibration-free table or a vibration-reducing steel plate. For preventive measures, contact your Nikon representative.
- Always use a power source that is free from electric noise or sudden voltage fluctuations. An improper power source may cause degradation, malfunction, or failure of the laser unit or PC.
- To prevent a failure, do not obstruct the air flow by blocking the cooling fans of the AX controller, PC, or laser unit. When the system is installed close to a wall, provide space of at least 15 cm between the wall and each unit of the system.
- Install the product so that the power cord to the AX controller can be unplugged immediately in case of an emergency.

System diagram

General Room and Equipment Recommendations:

Room size:

3 x 3 meter minimum recommended for Ax and AxR (when set up in L shape).
The doorway size must accommodate installation of required anti-vibration table.

Vibration:

5um pp @5Hz; 10um pp @10Hz, 10um pp @20Hz.
Steps should be taken to avoid high air currents. High volume/Low flow air handling is suggested.
Air conditioning vents directly above instrumentation should be avoided.

Internet:

An internet connection for the host PC is recommended for remote diagnostics, training, and/or support purposes.

OS accessibility:

Nikon will require administrative privileges to the Windows OS on Ax workstations. This is required for setup and diagnostic purposes and can be supervised by IT personnel, if necessary.

Other:

Room should be clean and dust free.
Dark wall, Blackout curtains and/or dimmable overhead lights are suggested:
A dark environment is recommended for confocal imaging.
If wireless phone service is not possible, a working telephone is strongly suggested in the equipment room.

Customer's Responsibilities for Site Preparation:

PC Table:

A separate table or cart is recommended for placement of the PC computer and monitor. This furniture should be large enough accommodate a PC keyboard, mouse and monitor. Nowadays often curved monitors are used. They are quite wide (like the HP Z38c is 90cm wide). Nikon strongly recommends never putting the PC on the active vibration isolation table.

Active Vibration Isolation Table:

A minimum sized 80 x 60 cm anti vibration table with Active Vibration dampening is required for optimal system performance.

Often a compressed air or nitrogen supply regulated in the range of 60-80 psi is required to accommodate the anti-vibration table.

Please note that a regulator supplied by the user is required when using nitrogen tanks.

Nikon requires active vibration isolation tables to meet performance standards of the instrument.

Tables not purchased from Nikon Instruments need to meet minimum requirements for size and vibration dampening. (Detailed later in this document)

Through holes on the tabletop are recommended in either English or Metric arrangements for Ax systems.

If not purchased through Nikon Instruments, the buyer is responsible if equipment is damaged or equipment does not meet performance standards due to incorrect mounting or vibration dampening.

An anti-vibration table will be delivered in a large wooden crate and can be quite heavy.

Make sure the building entrance and elevator are sufficient sized and rated to transport the table in the crate to the room.

Please note that no fan driven controllers, mouse, keyboard or monitor may be installed on the anti-vibration table, so take measures to arrange where they can be placed. Like a shelf on the wall for the fan driver controllers and a separate table or desk for the PC monitor, mouse and keyboard.

Although there are no official vibration tolerance specifications listed, we recommend an anti-vibration table in specific usage conditions.

It mainly depends on what kind of magnification the end user will use to image the sample.

Worst case is 100x objective lens with a super high magnification in AX scan (like 500x ~ 1000x).

In this case, we are looking for Super Resolution vibration tolerances.

This SR specs is listed as "the floor vibration is 50 um/sec or below (Class A level)"

Best case is use of 20x with 1x AX scan.

In this case, we actually even do not need an anti-vibration table as long as the floor is stable.

In practice the recommendation will be between 100 and 50 um/sec.

Criterion Curve (see Figure 1)	Max Level (1) micrometers/ sec,rms	Detail Size (2) microns	Description of Use
Workshop (ISO)	800	N/A	Distinctly feelable vibration. Appropriate to workshops and nonsensitive areas.
Office (ISO)	400	N/A	Feelable vibration. Appropriate to offices and nonsensitive areas.
Residential Day (ISO)	200	75	Barely feelable vibration. Appropriate to sleep areas in most instances. Probably adequate for computer equipment, probe test equipment and low-power (to 20X) microscopes.
Op. Theatre (ISO)	100	25	Vibration not feelable. Suitable for sensitive sleep areas. Suitable in most instances for microscopes to 100X and for other equipment of low sensitivity.
VC-A	50	8	Adequate in most instances for optical microscopes to 400X, microbalances, optical balances, proximity and projection aligners, etc.
VC-B	25	3	An appropriate standard for optical microscopes to 1000X, inspection and lithography equipment (including steppers) to 3 micron line widths.
VC-C	12.5	1	A good standard for most lithography and inspection equipment to 1 micron detail size.
VC-D	6	0.3	Suitable in most instances for the most demanding equipment including electron microscopes (TEMs and SEMs) and E-Beam systems, operating to the limits of their capability.
VC-E	3	0.1	A difficult criterion to achieve in most instances. Assumed to be adequate for the most demanding of sensitive systems including long path, laser-based, small target systems and other systems requiring extraordinary dynamic stability.

Notes:

(1) As measured in one-third octave bands of frequency over the frequency range 8 to 100 Hz.

(2) The detail size refers to the line widths for microelectronics fabrication, the particle (cell) size for medical and pharmaceutical research, etc. The values given take into account the observation that the vibration requirements of many items depend upon the detail size of the process.

Room Air Conditioning:

Air conditioning units must be designed to accommodate the heat load of the instrumentation (details provided in this document below).

Airflow in instrument rooms should be minimized whenever possible.

Nikon recommends high volume/low flow air handling for this purpose.

Air conditioning vents should not be placed such that airflow is directed towards the instrumentation (for example, do not place an air vent directly above a confocal system).

Confocal systems placed near or under air-conditioning vents will experience vibrations visible in image data. Upright systems are especially susceptible to such influences.

Chairs/Stools:

Nikon recommends adjustable chairs/stools that can be extended to approximately 80 cm high for use with microscopes on anti-vibration tables.

Accessories Table or Shelving:

A table large enough to accommodate accessory controllers or devices is recommended.

This is generally 2-5 units per system.

Placement of peripheral controllers on vibration isolation tables is not recommended.

Typical* accessories (sizes (WxHxD) in mm) include:

Microscope controllers (210 x 140 x 305)

Peripheral device controllers (300 x 100 x 200)

Epifluorescence Illuminators (125 x 163 x 263)

*Sizes are provided for guidance and will vary by configuration.

Alternatively, wall shelving around or near the instrumentation (maximum distance 3 meter, recommended distance 1meter, due to cord lengths) can be used for accessory equipment provided:

- The shelving is at least 35cm in depth.
- 2,5 cm of space is available behind the shelving or drilled holes are placed in the bottom of the shelving for power and controller cables to pass through.
- The shelves can accommodate at least 25kg of weight.

Please Note:

Ax system components are stackable, free standing on wheels and do not require shelving or tables for placement.

These components include controller, detector units, and laser units.

Power Requirements:

Room power requirements are the responsibility of the customer and should be met prior to installation date of the equipment.

Details on power requirements for the Ax systems are outlined in this document.

The AX system, which includes the microscope, PC, laser unit, and all accessories requires at minimum two 240V 16Amp outlets on separate circuits.

A minimum of total 6 receptacles wall outlets are suggested.

These outlets are ideally placed on the wall behind the location of the anti-vibration table and/or the PC furniture.

Nikon recommends plugging the laser unit into one circuit and all other devices into the second.

Input Voltage Specifics:

+/-5%; Frequency 50/60Hz +/- 5%; Harmonic Distortion: <3%. Power free of surges and harmonics.


Line Conditioning & Battery Backup:

Nikon strongly suggests using in-line conditioners with AX equipment.

Optionally, you may want to consider battery backup to protect the system in the event of a power failure, particularly to prevent data loss with long-term experiments.

Dimensions and Weights:

LUNS4-S6 Laser Unit	External Dimensions	Weight	400mm (W) x 210mm (H) x 530mm (D) 20kg (Lasers Included)
LUNF Laser Unit	External Dimensions	Weight	350mm (W) x 195mm (H) x 450mm (D) 25kg (Lasers Included)
AX/AXR Scan Head	External Dimensions	Weight	264mm (W) x 170mm (H) x 330mm (D) 12kg
AX Controller	External Dimensions	Weight	400mm (W) x 390mm (H) x 530mm (D) 28kg
AX Interlock Box	External Dimensions	Weight	210mm (W) x 50mm (H) x 142mm (D) 1kg
DUX-ST detector	External Dimensions	Weight	400mm (W) x 269mm (H) x 530mm (D) 22kg (Max)
DUX-VB Detector (5)	External Dimensions	Weight	400mm (W) x 269mm (H) x 530mm (D) 25kg (Max)

System Heat Loads and Ventilation Requirements

Component	Watt (Max)	BTU/min	BTU/hour
Ax Scan Head, Controller, DUX, Microscope	800	45.5	2730
Laser unit LUNS4 or S6	200	12	683
Laser unit LUN4	200	12	683
Host PC (Max)	1200	68.25	4095
Host Monitor	200	12	683

Laser Specifications, Classification for Ax Systems:

Laser Type:

Continuous Wave.

Laser Class:

Class 3B for LUN4/LUNS4/LUNS6/LUNF (<500mW)

Operating Wavelengths*: 405, 440, 457, 477, 488, 514, 532, 543, 561, 640, 647, 730.

*Wavelengths on a system will depend upon the laser configuration chosen; not all wavelengths are applicable on all systems.

Beam Diameter:

0.8mm.

Maximum Expected Exposure:

No direct exposure to coherent beam; beam is introduced via closed, interlocked optical path

Safety:

Laser heads are housed in a protective, light-tight cover.

For LUN4/S4/S6, there are mechanical safety shutters for each line and an AOTF that must be energized for laser light to reach the microscope.

For LUNF, redundant electrical circuits are interlocked and directly cut off power to each diode.

Laser is delivered to the scan head via a shielded optical fiber.

The microscope is mechanically interlocked via an interlock switch with a confocal controller.

A 100% reflection laser-safe beam splitter is used to direct light into the confocal optical path.

This laser product is designed and manufactured in compliance with the "Performance Standards for Light-Emitting Products" stipulated by the FDA in the U.S.A. and "Safety of Laser Products" stipulated by the IEC. Compliance with FDA standards is in accordance with Laser Notice No. 56. Users are requested to take all appropriate safety measures specified by the preceding standards, in accordance with local laws and regulations.

Use of this laser product must be supervised by a laser safety officer (LSO) and users must follow the LSO's instructions.

IEC Class 3B Laser Product

Wavelength: 400 to 700 nm
(continuous wave (CW) visible laser)

Beam divergence: 0.26 to 2.74 rad

Maximum power or energy output:
490 mW (CW visible laser)

Complies with FDA performance standards for laser products except for conformance with IEC 60625-1 Ed. 3., as described in Laser Notice No. 56, dated May 8, 2019.

IEC Class 4 Laser Product

Wavelength: 400 to 700 nm
(continuous wave (CW) visible laser)

Beam divergence: 0.26 to 2.74 rad

Maximum power or energy output:
600 mW (CW visible laser)



Complies with FDA performance standards for laser products except for conformance with IEC 60625-1 Ed. 3., as described in Laser Notice No. 56, dated May 8, 2019.


Suggested Room Layout for Ax Systems

Layout of the AX System (With the Laser Unit LUA-S / LUD-S Series):



Typical layout for right handed end users:

Layout in "L" Arrangement for smallest footprint:

Layout of the AX System (With the Laser Unit LUD-H Series):

Cautions

1. Setup of the system

This system shall be set up by trained Nikon service personnel. Assigning this task to any other personnel may result in improper setup, impairing system performance, and exposing eyes or skin to the laser beam, resulting in injury.

2. Prohibition of disassembly

Disassembling this system may result in electric shock or other accidents.

Do not attempt to disassemble any unit that constitutes this system.

In particular, disassembling any unit while the laser is active may cause laser light leakage.

3. Power supply ratings and power cords

Before plugging the power cords into the power outlets, check the power ratings to make sure that the input voltage and current capacity requirements of the product are met. Use of a nonconforming power source may result in unit malfunction, failure, or fire. Always use power cords that satisfy local safety regulations. Always use three-prong grounded power outlets.

Turn off the power switches of all the units before plugging or unplugging the power cords. To prevent breakage of electrical wires when plugging or unplugging power cords, do not pull it by the jacketed-cord part. Firmly grasp the plug instead.

4. Cable connections

All the cables in this system shall be connected by trained Nikon service personnel. However, you can connect the optional AX external trigger cable yourself.

- Do not disconnect or reconnect the cables that have been connected by Nikon service personnel. If any abnormal state is found on the cable connections, do not touch the cable and contact your nearest Nikon representative.
- Before connecting the AX external trigger cable, turn off the power to all the units. For turning off the power to the laser unit, see the instruction manual supplied with the laser unit.
- When plugging or unplugging the AX external trigger cable, do not pull it by the jacketed-cable part. Firmly grasp the plug instead; otherwise, the electrical wires may break.

5. Countermeasures for problems

If the system is not working properly, immediately terminate the software on the PC, turn off the power switches of all the units, unplug all the power cords from the power outlets, and then contact your nearest Nikon representative.

6. Avoid contact with water

If water comes in contact with a unit of the system, a short circuit might occur and cause malfunction. If water splashes on a unit, immediately turn off the power switch of the unit, and then wipe off moisture with a dry cloth or other suitable material. If water enters the unit, immediately stop using the system and unplug the power cords from the power outlets, and then contact your nearest Nikon representative.

7. Cautions against placement of units on the other units and cautions about relocation

Setting up of this system entails specific placement of the units on the other units. These units and the placement order are predetermined.

Because this system has been set up by trained Nikon service personnel, do not additionally

place any other unit or change the placement of the units.

If the units placed on the other units must be relocated, push the AX controller (equipped with wheels), which is the unit placed at the bottom. Carelessly pushing any of the units placed on the other unit(s) may cause it to drop, which may cause injury, or deformation, damage, or failure of the unit.

When moving the unit placed on the other unit(s) individually, always hold and lift the unit with both hands. When lifting the detector unit, in particular, hold on the bottom surface of the detector unit tightly; because, it is heavy.


Standards

<p>EU</p>	<p>CE marking</p> <ul style="list-style-type: none"> • Low Voltage Directive EN61010-1 (Safety requirements for electrical equipment for measurement, control and laboratory use) EN60825-1 (Laser safety) • EMC Directive EN61326-1 (EMC for electrical equipment for measurement, control and laboratory use) EN55011 Class A (EMI) EN61000-3-2 (Harmonics) EN61000-3-3 (Flicker) EN61000-4-2 (Electro-static discharge) EN61000-4-3 (Radiated immunity) EN61000-4-4 (Electrical fast transients) EN61000-4-5 (Surge) EN61000-4-6 (Conducted immunity) EN61000-4-8 (Power frequency magnetic) EN61000-4-11 (Dips and sags) 
<p>USA</p>	<p>FDA</p> <ul style="list-style-type: none"> • Notice 56 is applied. • FDA CDRH 21CFR Part1040.10 (Laser products performance) <p>FCC</p> <ul style="list-style-type: none"> • Subpart 15B ClassA (Unintentional Radiators) <p>Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.</p>
<p>Others</p>	<p>Canada ICES-003(A) / NMB-003(A) Australia AS/NZS CISPR11</p>

END.


FLIM adds an extra dimension to AX



Why is FLIM essential?

- Additional imaging parameter to gain more contrast
- Multiplexing
- Highly quantitative imaging method → FLIM-FRET
- QUAREP-LiMi has established a FLIM working group in 2024
 - Shows the demand and the potential for FLIM in confocal microscopy


Bellis - AX and FLIM images with Lambda D 20x objective lens Galvo Nyquist 1024 x 1024

FLIM applications

Identify fluorophores


Micro wellplate with **DASPI** (0.3 ns), **Coumarin6** (2.1 ns), **Fluorescein** (4.1 ns), **Atto645** (5.3 ns)

Remove background
Using timegating


2D Gel with different Cy2 labeled protein spots. In the right part the fast scattering and long background lifetime components are removed.

Autofluorescence imaging
through enhanced contrast


Long lifetime autofluorescence (green) from tumor and stromal cells versus fast fluorescence from collagen (blue).

Tissue differentiation by
fluorophore separation


Lung tissue cell types separated based on the lifetime

FRET to study interactions

GFP lifetime decrease in cytoplasm vesicles due to FRET with mRFP.

Material science

Single nanodiamonds with fluorescent NV-defect sites.

NEW state of the art laser technology

Former laser head

- Large footprint
- Elliptical beam shape 3-5mm
- Dual mode (pulsed and cw)
- No beam blanking
- Flexible repetition rate

New laser head - Launched in 2023

- Small cubic footprint (like OBIS laser)
- Round beam shape 1mm in diameter
- Dual mode (pulsed and cw)
- Beam blanking
- Flexible repetition rate

Shared Laser Unit (SLU) for Nikon AX

NEW!

Cost effective

Neat solution

Upgradability

Usability

Shared Laser Unit (SLU) for Nikon AX

Cost effective

One laser bed for
AX and FLIM
modalities

Neat solution

streamline
implementation

Upgradability

Configuration in
Building blocks
style

Usability

Lasers can be used for
each AX detectors

SLU: What's inside the box


- Switchable mode laser unit: Continuous Wave (CW) and pulsed for Confocal and FLIM respectively
- Dichroic displacer mirror for easy alignment
- AOTF guaranteeing fast switch between laser line and laser power stability and linearity
- Guaranteed laser power 15mW
- One Obis laser 561nm (not switchable to pulse)
- Flexible repetition rate
- On-field maintenance


Hardware is controlled by Nikon and unit is neatly integrated into Nikon
Software: NIS-element

Neat solution – Seamless integration in NIS

Neat solution – Seamless integration in NIS

The screenshot displays the Nikon NIS-Elements AR software interface, showcasing a seamless integration of various microscopy components. The interface is divided into several key sections:

- Left Panel (ROI Statistics):** Displays a histogram and a line graph showing the intensity profile of a selected region of interest (ROI). Below the graphs, a table provides statistical data for the ROI.
- Center-Left Panel (ROI Statistics Table):**

Property	Value	Unit
Area [μm²]	8762.58	N/A
Mean First Moment	1.04	N/A
Min First Moment	0	N/A
Max First Moment	18.77	N/A
Sum First Moment	189218.24	N/A
SD of First Moment	1.67	N/A
Subst. Background	N/A	N/A
- Center-Right Panel (Acquisition):** Shows the acquisition settings for the current experiment, including the camera (AX-FIM A FCS), channels (GFP, TRIC2), and various parameters like dwell time (0.4 sec) and frame rate (0.506 fps).
- Right Panel (Microscopy):** Displays the microscope's configuration, including the objective lens (1.8 NA, 400 mm), field of view (177 μm x 177 μm), and other technical specifications.
- Main View:** Shows a live microscopy image of a cell culture, with a green ROI overlaid on a purple-stained cell. The image is displayed in a large window, allowing for detailed observation and manipulation.

SLU configuration – LU PQ 4/6

4 lines 405/488/561/640

6 lines 405/445/488/515/561/640

Upgrade 594nm instead of 561 possible
LU PQ 4 Lines to 6 Lines on field

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EU Declaration of Conformity

This European Declaration of Conformity is issued under the sole responsibility of the manufacturer.

MANUFACTURER	
Name of Company	Address
Nikon Corporation	471, Nagaodai-cho, Sakae-ku, Yokohama, Kanagawa 244-8533 Japan

AUTHORIZED REPRESENTATIVE		
Name of Company	Address	Telephone
Nikon Europe B.V.	Stroombaan 14, 1181 VX Amstelveen, The Netherlands	+31-20-7099-000

PRODUCT IDENTIFICATION	
Product / Trade Name	Product Code
AX-SHS, AX-SHR	MHA52000, MHA52100

STANDARDS / COMMON SPECIFICATIONS	
Standards	Common Specifications
EN 61010-1:2010/A1:2019/AC:2019-04 EN 60825-1:2014+A11:2021 EN 61326-1:2013 Class A EN IEC 63000:2018	N/A

Nikon Corporation declares that the above-mentioned products meet the provision of the following EU legislation:

- Low Voltage Directive (2014/35/EU)
- EMC Directive (2014/30/EU)
- RoHS Directive (2011/65/EU, 2015/863/EU)
- Regulation (EU) 2023/1542

COMPANY REPRESENTATIVE: HARUYUKI YAMASHITA

TITLE: Department Manager
Quality Assurance Department

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